


Arthroscopy




MAHE International, Inc.

Made in Germany

CE / ISO 9001

Arthroscopes and Accessories



Small Joint Arthroscopes

Ø 2.7 mm, working length 145 mm

	Non-Autoclavable	Autoclavable
0°	50.002.00	50.002.00A
30°	50.002.30	50.002.30A
70°	50.002.70	50.002.70A

AUTOCLAVE
134° C / 273° F

Arthroscopes


Ø 2.7 mm, working length 176 mm

	Non-Autoclavable	Autoclavable
0°	50.003.00	50.003.00A
30°	50.003.30	50.003.30A
70°	50.003.70	50.003.70A

Arthroscopes

Ø 4 mm, working length 176 mm

	Non-Autoclavable	Autoclavable
0°	50.004.00	50.004.00A
30°	50.004.30	50.004.30A
70°	50.004.70	50.004.70A



Sheaths and Trocars

Ø 2.7 mm, working length 145 mm

Sheath with 1 fixed stopcock	50.021.02
Sheath with 2 fixed stopcocks	50.022.02
Sheath with 1 rotatable stopcock	50.031.02
Sheath with 2 rotatable stopcocks	50.032.02
Trocar, pointed	50.035.02
Trocar, blunt	50.036.02

Ø 2.7 mm, working length 176 mm

Sheath with 1 fixed stopcock	50.021.03
Sheath with 2 fixed stopcocks	50.022.03
Sheath with 1 rotatable stopcock	50.031.03
Sheath with 2 rotatable stopcocks	50.032.03
Trocar, pointed	50.035.03
Trocar, blunt	50.036.03


Ø 4 mm, working length 176 mm

Sheath with 1 fixed stopcock	50.021.04
Sheath with 2 fixed stopcocks	50.022.04
Sheath with 1 rotatable stopcock	50.031.04
Sheath with 2 rotatable stopcocks	50.032.04
Trocar, pointed	50.035.04
Trocar, blunt	50.036.04

Ø 5.5 mm, working length 176 mm

Sheath with 1 fixed stopcock	50.021.05
Sheath with 2 fixed stopcocks	50.022.05
Sheath with 1 rotatable stopcock	50.031.05
Sheath with 2 rotatable stopcocks	50.032.05
Trocar, pointed	50.035.05
Trocar, blunt	50.036.05

Hand Instruments with knurled handle - Ø 4 mm




Probe, graduated
3.5 mm
50.050.01



Scalpel Knife
50.050.02



Banana Knife,
double-edged
50.050.03



Banana Knife,
tooth-edged
50.050.04



Menisctome
3 mm
50.050.05



Menisctome
5 mm
50.050.06



Hook Knife
50.050.07



Retrograde Knife
50.050.08



Ring Curette, 3mm,
15° angled
50.050.09



Rosette Knife,
tooth-edged
50.050.10



Curette, 3 mm,
15° angled
50.050.11

MAHE International, Inc. 468 Craighead Street Nashville, TN 37204 USA

Website: MahelInternational.com


Toll Free (800) 294-7946 Phone (615) 269-7256 Fax (615) 269-4605 Email info@MahelInternational.com

ART-01

1/2005



Hand Instruments with adaptable or triangular handle




Probe, graduated

50.051.01 adaptable
50.052.01 triangular



Probe, graduated

50.051.02 adaptable
50.052.02 triangular



Curette, 3mm,

50.051.03 adaptable
50.052.03 triangular



Ring Curette, 3 mm,

50.051.04 adaptable
50.052.04 triangular



Ring Curette, 5 mm,

50.051.05 adaptable
50.052.05 triangular



Curette

50.051.06 adaptable
50.052.06 triangular



Cartilage Files

50.051.07 adaptable
50.052.07 triangular



Cartilage Files

50.051.08 adaptable
50.052.08 triangular



Meniscotome

50.051.09 3 mm adaptable
50.052.09 3 mm triangular
50.051.10 5 mm adaptable
50.052.10 5 mm triangular



Smillie Meniscotome

50.051.11 adaptable
50.052.11 triangular



Rosette Knife, 3 mm,

50.051.12 adaptable
50.052.12 triangular



Rosette Knife, 5 mm,

50.051.13 adaptable
50.052.13 triangular



Scalpel Knife

50.051.14 adaptable
50.052.14 triangular




Bayonet Knife

50.051.15 adaptable
50.052.15 triangular



Banana Knife,

50.051.16 adaptable
50.052.16 triangular



Banana Knife,

50.051.17 adaptable
50.052.17 triangular



Sickle Knife

50.051.18 adaptable
50.052.18 triangular



Hook Knife

50.051.19 adaptable
50.052.19 triangular



Hook Knife

50.051.20 adaptable
50.052.20 triangular



Retrograde Knife

50.051.21 3 mm adaptable
50.052.21 3 mm triangular

Adaptable
Handle only
50.051.00

MAHE International, Inc. 468 Craighead Street Nashville, TN 37204 USA


Website: MahelInternational.com















ART-02

1/2005

Toll Free (800) 294-7946 Phone (615) 269-7256 Fax (615) 269-4605 Email info@MahelInternational.com

Scissors and Punches



	 Hook Scissors	 Hook Punch with scoop	 Hook Punch, blunt
 Straight	Ø 2,7 mm 51.000.04 3,4 mm 51.000.06	Ø 2,0 mm 51.001.02 2,7 mm 51.001.04 3,4 mm 51.001.06	Ø 2,0 mm 51.002.02 2,7 mm 51.002.04 3,4 mm 51.002.06
 30° right	Ø 2,7 mm 51.000.12 3,4 mm 51.000.14	Ø 2,7 mm 51.001.12 3,4 mm 51.001.14	Ø 2,7 mm 51.002.12 3,4 mm 51.002.14
 30° left	Ø 2,7 mm 51.000.20 3,4 mm 51.000.22	Ø 2,7 mm 51.001.20 3,4 mm 51.001.22	Ø 2,7 mm 51.002.20 3,4 mm 51.002.22
 45° right	Ø 2,7 mm 51.000.28 3,4 mm 51.000.30	Ø 2,7 mm 51.001.28 3,4 mm 51.001.30	Ø 2,7 mm 51.002.28 3,4 mm 51.002.30
 45° left	Ø 2,7 mm 51.000.36 3,4 mm 51.000.38	Ø 2,7 mm 51.001.36 3,4 mm 51.001.38	Ø 2,7 mm 51.002.36 3,4 mm 51.002.38
 7° upward	Ø 2,7 mm 51.000.44 3,4 mm 51.000.46	Ø 2,7 mm 51.001.44 3,4 mm 51.001.46	Ø 2,7 mm 51.002.44 3,4 mm 51.002.46
 7° downward	Ø 2,7 mm 51.000.52 3,4 mm 51.000.54	Ø 2,7 mm 51.001.52 3,4 mm 51.001.54	Ø 2,7 mm 51.002.52 3,4 mm 51.002.54
 15° upward	Ø 2,7 mm 51.000.60 3,4 mm 51.000.62	Ø 2,7 mm 51.001.60 3,4 mm 51.001.62	Ø 2,7 mm 51.002.60 3,4 mm 51.002.62
 15° downward	Ø 2,7 mm 51.000.68 3,4 mm 51.000.70	Ø 2,7 mm 51.001.68 3,4 mm 51.001.70	Ø 2,7 mm 51.002.68 3,4 mm 51.002.70
 7° jaw upward	Ø 2,7 mm 51.000.76 3,4 mm 51.000.78	Ø 2,7 mm 51.001.76 3,4 mm 51.001.78	Ø 2,7 mm 51.002.76 3,4 mm 51.002.78
 15° jaw upward	Ø 2,7 mm 51.000.84 3,4 mm 51.000.86	Ø 2,7 mm 51.001.84 3,4 mm 51.001.86	Ø 2,7 mm 51.002.84 3,4 mm 51.002.86

MAHE International, Inc. 468 Craighead Street Nashville, TN 37204 USA

Website: MahelInternational.com









Toll Free (800) 294-7946 Phone (615) 269-7256 Fax (615) 269-4605 Email info@MahelInternational.com


ART-03

1/2005



Scissors, Punches, Side Biter and Accessories

Ø 3.4 mm	 MAHE Wide Jaw Punch, large	 MAHE Wide Jaw Punch, medium
 Straight	<u>51.050.00</u>	<u>51.051.00</u>
 30° right	<u>51.050.02</u>	<u>51.051.02</u>
 45° right	<u>51.050.03</u>	<u>51.051.03</u>
 30° left	<u>51.050.04</u>	<u>51.051.04</u>
 45° left	<u>51.050.05</u>	<u>51.051.05</u>
 15° upward	<u>51.050.10</u>	<u>51.051.10</u>



Hook Scissors 45°

3.4 mm right	<u>51.008.02</u>
3.4 mm left	<u>51.008.04</u>




Hook Punch with scoop 45°

3.4 mm right	<u>51.009.02</u>
3.4 mm left	<u>51.009.04</u>




90° Wide Jaw Punch

3.4 mm right	<u>51.053.12</u>
3.4 mm left	<u>51.053.14</u>




Retrograde Wide Jaw Punch

3.4 mm right	<u>51.055.02</u>
3.4 mm left	<u>51.055.01</u>




Grasping Forceps, without ratchet

2.7 mm	<u>51.060.00</u>
3.4 mm	<u>51.060.02</u>



Grasping Forceps 1 x 2 Teeth, without ratchet

2.7 mm	<u>51.060.10</u>
3.4 mm	<u>51.060.12</u>



Grasping and Biopsy Forceps, without ratchet

2.7 mm	<u>51.060.20</u>
3.4 mm	<u>51.060.22</u>

Grasping Forceps, with ratchet


2.7 mm	<u>51.060.04</u>
3.4 mm	<u>51.060.06</u>

Grasping Forceps 1 x 2 Teeth, with ratchet


2.7 mm	<u>51.060.14</u>
3.4 mm	<u>51.060.16</u>

Grasping and Biopsy Forceps, with ratchet

2.7 mm	<u>51.060.24</u>
3.4 mm	<u>51.060.26</u>



51.080.00



	Ø 3.2 mm	Ø 4.5 mm
Drainage Cannula with 1 tap	<u>52.150.32</u>	<u>52.150.45</u>
Trocar, oblique tip	<u>52.151.32</u>	<u>52.151.45</u>
Obturator, blunt	<u>52.152.32</u>	<u>52.152.45</u>

MAHE International, Inc. 468 Craighead Street Nashville, TN 37204 USA

Website: MahelInternational.com

ART-04

1/2005

Toll Free (800) 294-7946 Phone (615) 269-7256 Fax (615) 269-4605 Email info@MahelInternational.com