


180 W Xenon light source



- Electronically controlled light intensity
- Audible overheating alarm
- Easy lamp replacement


10.912.00

Light Source Specifications


Lamp Type	OSRAM XBOR R180/45
Lamp Power	180 W
Lamp Color Temperature	6000° K
Lamp Life	500 hours minimum
Weight	5.5 kg
Dimensions	295 (W) x 135 (H) x 253(D)

Universal Light cable

Ø	Working length	
3.5 mm	1800 mm	11.035.18
3.5 mm	2300 mm	11.035.23
3.5 mm	3000 mm	11.035.30
4.8 mm	1800 mm	11.048.18
4.8 mm	2300 mm	11.048.23
4.8 mm	3000 mm	11.048.30




Light Cable - Light Source - Adapter



Storz		11.001.00
Olympus		11.001.02
Wolf / Dyonics		11.001.04
ACMI / British Standard		11.001.06
Pentax		11.001.08
Fuji		11.001.10

Light Cable - Endoscope - Adapter



Storz / Olympus		11.001.20
Wolf		11.001.22
ACMI / British Standard		11.001.24

MAHE International, Inc. 468 Craighead Street Nashville, TN 37204 USA

Website: MahelInternational.com

Toll Free (800) 294-7946 Phone (615) 269-7256 Fax (615) 269-4605 Email info@MahelInternational.com